bewell[®] connect




MyOxy

CE
0197

BW-OX1

By **Visiomed**

SN Serial number

SN:

Year Month Day Serial number

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MyOxy

CE
0197

THIS PRODUCT RESPECTS THE REQUIREMENTS OF FCC, ROHS, IEC60601-1, IEC60601-1-2, IEC60601-1-11, ISO80601-2-61, ISO10993-1, ISO10993-5 AND ISO10993-10 STANDARDS.

CONTENTS

1. WARNING	16
2. INFORMATIONS	17
2-1. OXYGEN SATURATION (SpO ₂)	17
2-2. PULSE	18
3. FEATURES	19
4. USING THE OXIMETER.....	20
4-1. DESCRIPTION	20
4-2. PRECAUTIONS BEFORE USE	20
4-3. RECOMMENDATIONS BEFORE TAKING MEASUREMENTS	21
4-4. DOWNLOADING THE APPLICATION	21
4-5. HOW TO TAKE MEASUREMENTS	21
5. CLEANING AND MAINTENANCE	22
6. TECHNICAL SPECIFICATIONS	22
7. TROUBLESHOOTING	23

Find out how to install and use the BewellConnect® application by visiting our website www.bewell-connect.com/install

The manufacturer reserves the right to modify the technical characteristics of the product without prior notice.

Dear Customer,

Thank you for purchasing this BewellConnect® digital pulse oximeter. We hope you will be fully satisfied with it and we recommend reading this user's guide carefully so that you can measure your arterial oxygen saturation (SpO₂) and your pulse rate (PRbpm) effectively.

1. WARNING

- Only use this device for the purpose for which it was designed as outlined in this guide.
- This device is adapted for personal use at home.

- Use this device in a temperature range between 5 and 40°C.
- Do not expose this device to extreme temperature conditions (over 70°C or minus 25°C).
- Do not use this device in over 93% of relative humidity.
- This device must always be kept in a clean and dry place.
- Do not expose the device to sunlight or water.
- Do not expose the device to electric shocks.
- Do not drop the device.
- Follow the maintenance instructions given in this guide.
- Do not try to open the device. In the event of malfunctions, contact your retailer.
- Keep out of reach of children.
- The measurement results are given as a guide. If you have any doubts about your results, please contact your doctor.
- Stop using the device in the event of an anomaly or malfunction.
- This device is not designed for use by people (including children) with reduced physical, sensory or mental capacities, or by people lacking experience or knowledge, unless they are assisted by a person who is responsible for their safety and supervision or who has been given instructions about using this device.

2. INFORMATIONS

2-1. OXYGEN SATURATION (SpO₂)

Oxygen saturation measures the percentage of hemoglobin binding sites in the bloodstream that is occupied by oxygen. Many factors can modify hemoglobin's affinity for oxygen and reduce oxygen saturation in the blood. Measuring oxygen saturation is vital for individuals who may be at risk of desaturation, a condition that could result in the need for medical treatment to address the deficiency. A finger pulse oximeter can help monitor oxygen saturation.

ANALYSIS OF RESULTS

Warning: This table DOES NOT apply to individuals already suffering from certain illnesses (e.g. asthma, heart disease, respiratory diseases) or who are at altitudes above 1500 metres. If you are already suffering from a disease, always consult your doctor to analyze your results. If you have any doubts about your results, contact your doctor.

SpO ₂ (peripheral capillary oxygen saturation) result as a %	Diagnosis
99 - 94	Normal level
94 - 90	Reduced level (medical visit recommended)
< 90	Critical level : consult a doctor immediately or call for emergency help

FALL IN PERIPHERAL CAPILLARY OXYGEN SATURATION DUE TO ALTITUDE

Please note: The following table gives information about the effects of different altitudes on oxygen saturation levels as well as their effects on the human body. The following table DOES NOT apply to individuals already suffering from certain illnesses (e.g. asthma, heart disease, respiratory diseases). In people already suffering from illness, the symptoms (e.g. hypoxia) may already occur at low altitudes.

Altitude	SpO ₂ (peripheral capillary oxygen saturation) result as a %	Consequences for the individual
1500 - 2500 m	> 90	No mountain sickness (in general)
2500 - 3500 m	~ 90	Mountain sickness, adaptation recommended
3500 - 5800 m	< 90	High probability of mountain sickness, adaptation imperative
5800 - 7500 m	< 80	Severe hypoxia, only a stay for a limited period of time is possible
7500 - 8500 m	< 70	Immediate life-threatening danger

2-2. PULSE

Warning: the results below show the average for a person in good health whose measurements are taken in a calm environment.

An individual may have a higher or lower pulse rate as a result of numerous factors (illness, regular exercise, medical treatment, etc.).

In the event of irregularities in your results (lower or higher pulse rate than normal), we advise you to consult your doctor.

AVERAGE RESULTS

Average result for pulse rate (in beats per minute)	User range
120	Infant
100	Child aged 1 - 8
60 - 80	Adult

ASSESSMENT OF RESULTS

Normal pulse rate result (in beats per minute)	High pulse rate result (in beats per minute)	Low pulse rate result (in beats per minute)	Age range
120 - 150	> 150	< 120	Infant
80 - 150	> 150	< 80	Child aged 1 - 5
60 - 120	> 120	< 60	Child aged 5 - 12
60 - 105	> 105	< 60	Adolescent
60 - 80	> 80	< 60	Adult
90	-	-	Elderly person

3. FEATURES

The digital pulse oximeter is a portable device designed for occasional monitoring of arterial oxygen saturation (as a % of SpO₂) and pulse rate (PR) in beats per minute (bpm) on adult and child patients, in hospitals, health centres and at home. It should not be used for continuous monitoring (it does not include an alarm signalling an abnormal result), but it can provide continuous measurements for as long as the patient's finger remains in the device.

- Simple and practical to use.
- A compact product that is light and easy to carry around.
- LED display.

4. USING THE OXIMETER

4-1. DESCRIPTION

See *schema A* page 2

- 1 Readings screen
- 2 USB port
- 3 SpO₂ level
- 4 Pulse rate
- 5 Bluetooth LED
- 6 Signal indicator
- 7 Battery indicator

The Bluetooth LED flashes red when searching the Bluetooth signal and is fixed red when the connection is made.

If the screen shows "??", it means the signal is unstable. You need to take the reading again.

The battery indicator flashes when the battery is low. The device must be recharged using the USB cable (See *schema B* page 2).

4-2. PRECAUTIONS BEFORE USE

- Do not use this device near an Electrosurgical Unit (ESU), since it may function incorrectly as a result.
- The pulse oximeter must be placed so that the pulse can be recorded correctly and the main objective is to evaluate the oxygen saturation (SpO₂) level. Check that nothing is interfering with the readings.
- Do not use the device in an MRI (Magnetic Resonance Imaging) or CT (Computed Tomography) environment.
- Do not use this device in situations requiring the presence of an alarm. The device is not fitted with an alarm and is not designed for continuous monitoring.
- Do not use the device in an explosive atmosphere.
- The device is designed to help evaluate the patient's health. It must be used in conjunction with other methods to evaluate symptoms and clinical signs.
- Make frequent checks (every 30 minutes) on the application zone of the device sensor (skin sensitivity, position).
- The device is not designed to be used while the patient is being transported outside the health centre.
- The device must not be used near or on top of another device.
- The device must not be used when it is charging.
- Do not insert anything into the device except the USB cable during charging.
- Respect the temperature ranges and humidity levels shown in the specifications.

Inaccurate readings may result from :

- Nail-varnish (colorimetric interference) or artificial nails
- High levels of dysfunctional hemoglobin (carboxyhemoglobin or methemoglobin)
- Intravascular dyes (indocyanine green or methylene blue)
- Bright light. If necessary, protect the sensor positioning area
- Excessive movement of the patient
- Defibrillators and interference from high-frequency electrosurgical devices
- Venous pulse
- Placing the sensor on a patient with a blood-pressure cuff, an arterial or intravascular catheter
- High blood pressure, severe vasoconstriction, acute anemia or hypothermia
- Heart failure or stress
- Poor pulse quality (low perfusion)
- Low level of hemoglobin.

4-3. RECOMMENDATIONS BEFORE TAKING MEASUREMENTS

- When using the device for the first time, make sure that the battery has been charged via the USB port. Connect the USB cable of the device to a computer to charge it. Once the device has been charged, you can disconnect the USB cable.
- Avoid eating, smoking or taking exercise for 30 minutes before measuring your oxygen saturation.
- Remain calm for 15 minutes before the measurement.
- Keep still and do not speak during the measurement.

4-4. DOWNLOADING THE APPLICATION

The following information may be updated.

Download the BewellConnect® application to your telephone or tablet :

- via the App Store or Google Play
- by scanning the QR Code on the side of the box.

Then click on the MyOxy icon.

The results are immediately transferred by the application via Bluetooth 4.0 Smart Ready.

4-5. HOW TO TAKE MEASUREMENTS

See *schema C* page 2

Insert your finger (generally the index finger) into the device. The fingernail must be facing upwards. The device turns itself on automatically when a finger is inserted. The readings appear on the device screen. When the result flashes, the data have been transferred to the application.

Important : the data will only be transferred to the application if the application is open while the device is being used and Bluetooth 4.0 Smart Ready is activated.

The device switches off automatically after eight seconds.

If it is difficult to obtain a reading, try with the other hand or with another finger on the same hand.

5. CLEANING AND MAINTENANCE

-To clean the silicone inside the oximeter that comes into contact with the finger, use 70° rubbing alcohol. Apply the alcohol using a soft clean cloth. Also use the rubbing alcohol to clean the finger being measured, before and after use.

- Do not pour or spray liquid on to the surface of the oximeter. Make sure that liquid does not come into contact with the device openings. Leave the oximeter to dry before using again.

- Do not sterilize the device in an autoclave, with ethylene oxide or by immersing it in liquid. The device must not be sterilized.

6. TECHNICAL SPECIFICATIONS

Product name	BewellConnect® MyOxy
Model	BW-OX1
Classification	Classe IIa Applied part type BF IPX22 Continuous use
Display type	LED screen
LED Specifications	Wavelength: Red = 660 ±3nm - Infrared = 905 ±10nm Radiant power : Red = 3.2mW - Infrared = 2.4mW
Dimensions	32 x 32 x 62 mm
Weight	40 g
Display and measuring ranges	SpO ₂ : Display range : 35 - 99% / Measuring range : 70 - 99% Pulse: Display and measuring ranges: 30 - 250 bpm
Accuracy	SpO ₂ : 70 - 100% : ± 3% Pulse : ± 2 bpm
Resolution	SpO ₂ : 1% Pulse : 1 bpm
Normal working conditions	Temperature : 5 - 40°C Humidity : 15-93%RH

Storage conditions	Temperature : -25 - 70°C Humidity : ≤93%RH Atmospheric pressure : 86-106Kpa
Technology	Bluetooth 4.0 Smart Ready
Battery	Lithium 280mAh
Autonomy	2 hours in continuous use or 400 separate readings
Automatic shut-down	8 seconds
Bluetooth transmission distance	10m
Frequency	2400 - 2483.5MHz
USB	USB port 2.0
Accessories provided	Pulse oximeter, instruction manual, quick-start guide, protective pouch, USB cable

7. TROUBLESHOOTING

If one of the following issues arises when using your MyOxy, consult this trouble-shooting guide. If the problem continues, please contact our customer service department.

ISSUE	POSSIBLE CAUSE	SOLUTION
SpO ₂ or PRbpm do not appear correctly	Your finger must be correctly inserted.	<i>Change the position of your finger and try again.</i>
	SpO ₂ reading too small to be measured.	<i>Try several times. If you have any doubts about a problem related to the product, seek medical attention to obtain an accurate diagnosis.</i>
SpO ₂ or PRbpm are unstable	Your finger may not be inserted deeply enough.	<i>Change the position of your finger and try again.</i>
	The patient's finger is trembling or the patient is moving.	<i>Try not to move.</i>
Impossible to turn on the device	Battery low	<i>Recharge the device with the USB cable</i>
	Oximeter is damaged	<i>Contact customer support</i>

bewell[®] connect

The indicator lights turn off	The device switches itself off automatically after eight seconds when a signal is not detected.	Normal
Error3	Low battery	Recharge the device with the USB cable
	Oximeter is damaged	Contact customer support
Error7	Low battery	Recharge the device with the USB cable
	Oximeter is damaged	Contact customer support
No data transferred to your telephone/tablet	Bluetooth connection deactivated	Activate the Bluetooth on your telephone/tablet
	Telephone or tablet incompatible	Check that your telephone or tablet has Bluetooth 4.0 Smart Ready

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EN: Visiomed® will repair or replace this product free of charge in the case of defective parts or manufacturing defects, in accordance with the conditions mentioned below as follows:

DURATION: 24 MONTHS RETURN TO WORKSHOP

LIMITS AND EXCLUSIONS: This guarantee concerns only the original final purchaser. A purchase invoice, or another proof of purchase, with this guarantee card will be required to obtain an after-sales service, in accordance with this guarantee. This guarantee card will not be extended to another person only the original final purchaser. This guarantee becomes void if the serial numbers on the product are modified, replaced, illegible, absent, or if repair has been carried out by a service not approved, including the user.

This guarantee covers only the defects of the material or parts, occurring during normal use of the product. It does not cover the damage caused during the transport of the apparatus, causes due to repairs being carried out by the distributor, by any modifications undertaken, any connection of equipment not approved by Visiomed®, or causes contrary to those written in the user manual or notice. Moreover, the present guarantee does not cover damage due to falls, bad handling, bad installations, damage by fire, floods, lightning, or any other natural disaster. This guarantee does not cover the packing of the material, the accessories, the defects caused by commercial exposure of the product, show room, sale space, demonstration etc... Normal maintenance, cleaning and the replacement of parts where wear is normal, are not covered by the terms of this guarantee. Visiomed® and its representatives and agents will not in any case be held responsible for any damage and consecutive damages due to the mishandling of this product. This guarantee is the only valid one at Visiomed®, any

other guarantee (commercial guarantee) except this one will not be taken into account.

IMPORTANT: During the guarantee period if you are dissatisfied with the repairs of this product, please contact the Visiomed® customer service.

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EN: *Disposal of Waste Equipment by Users in Private Households in the European Union.*

The symbol on the product or on its packaging indicates that this product must not be disposed of with your other household waste. Instead, it is your responsibility to dispose of your waste equipment by handing it over to a designated collection point for the recycling of waste electrical and electronic equipment. The separate collection and recycling of your waste equipment at the time of disposal will help to conserve natural resources and ensure that it is recycled in a manner that protects human health and the environment. For more information about where you can drop off your waste equipment for recycling, please contact your local city office or your household waste disposal service.

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EN: *BewellConnect® is a brand registered by Visiomed Group SA. The constructor keeps the right to modify the photographs and characteristics without previous notice. Registered design. The logos and trademarks are the property of their respective owners.*

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Made in P. R. C.
Nexmed Technology Co., LTD.
2nd Floor of N°.1 building, Jia An
Technological Industrial Park
67 District, Bao An
518101 Shenzhen - China

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Visiomed SAS France