bewell[®] connect

Distribué par / Distributed by
www.visiomed-lab.com

Marché Européen /
European Market
Visiomed SAS France
Service clients / Customers service
8, avenue Kléber
75116 Paris - France
Tel : +33 892 350 366 (0,34€/mn)
Fax : +33 1 44 17 93 10
contact@visiomed-lab.fr

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BW-TS1

By **Visiomed**

This **MyTens** device complies with all the applicable standards and regulations related to exposure to electromagnetic fields.

The device complies with all the standards relating to Class IIa electrical medical devices and to devices which use electrical stimulation for use at home.

Electromagnetic compatibility (EMC) : The **MyTens** device complies with safety standards EMC EN 60601-1-2. It is designed to be used at home or in hospitals.

CONTENTS

| | |
|--|-----------|
| 1. WARNINGS | 33 |
| 2. TRANSCUTANEOUS NEUROSTIMULATION INFORMATION | 35 |
| 3. CHARACTERISTICS | 36 |
| 4. USE | 36 |
| 4.1 DESCRIPTION | 36 |
| 4.2. PRECAUTIONS BEFORE USE | 36 |
| 4.3 PRESENTATION OF THE PROGRAMMES | 37 |
| 4.3.1. FORMS OF WAVES | 37 |
| 4.3.2. TENS PROGRAMMES | 38 |
| 4.3.3. ELECTRICAL MUSCLE STIMULATION (EMS) PROGRAMMES | 40 |
| 4.4. INDICATIONS BEFORE USE | 43 |
| 4.5. DOWNLOADING THE APPLICATION | 44 |
| 4.6. STARTING AND USING MYTENS | 44 |
| 5. STORAGE AND MAINTENANCE | 48 |
| 5.1. STORAGE OF THE ELECTRODES | 48 |
| 5.2. CLEANING | 48 |
| 6. TROUBLESHOOTING GUIDE | 49 |
| 7. AVOIDING MALFUNCTIONS | 49 |
| 8. TECHNICAL CHARACTERISTICS | 50 |
| 9. PLACING THE ELECTRODES | 51 |
| 9.1. POSITION OF THE ELECTRODES-TENS PROGRAMMES | 51 |
| 9.2. POSITION OF THE ELECTRODES-EMS PROGRAMMES | 53 |

Find out how to install and use the BewellConnect® app by visiting our website at www.bewell-connect.com/install

The manufacturer reserves the right to change the product's technical specifications without prior notice.

Dear customer,

Thank you for buying the BewellConnect® electro-stimulator. We hope you get the best use out of it, and we recommend that you read these instructions carefully so that you can use it as effectively as possible. To gain the best treatment from your **MyTens** device, please follow the treatment instructions in this user manual when using the device.

For specific information about vein, muscle or other problems, PLEASE CONTACT YOUR DOCTOR.

1. WARNINGS

- Please read this manual carefully and always follow the treatment instructions.
- Only use this device for its intended purpose as described in these instructions.
- This device may be used for personal use at home.
- Use this device in an ambient temperature range of 5 to 40°C.
- Do not expose this device to extreme temperature conditions > 50°C or < -10°C.
- Do not use this device at a relative humidity of over 85%.
- This device must always be placed in a clean and dry place.
- Do not expose this device to sunlight or water.
- Do not expose this device to electric shocks.
- Never drop the device.
- Follow the maintenance instructions specified in this manual.
- Do not attempt to open the device. In case of problems, contact your dealer.
- This is a medical device. Keep out of the reach of children.
- Discontinue use of the device in case of anomalies or malfunction.
- Electronic medical equipment requires special precautions regarding electromechanical compatibility. It must be installed and used in accord with EMC (ElectroMagnetic Compatibility) information.
- This device must not be used in an environment of strong electromagnetic interference: Near a TV, a microwave oven or a mobile phone in use, etc.
- This device is not designed to be used by persons (including children) whose physical, sensory or mental capabilities are reduced, or persons without experience or knowledge, unless they have been able to benefit, by the intermediary of a person responsible for their safety, from surveillance or prior instructions concerning use of the device. It is possible that they will not be able to use it in accordance with the instructions of this user manual and be disturbed by the treatment.
- This product is not designed for use on children.
- Do not modify the device or the electrodes. This could cause a malfunction.
- This device is designed for use by a single person and one person only.

Do not treat the the following parts or areas :

- Do not apply stimulation on your torso. Indeed, the introduction of electric current on this area can cause heart rhythm disturbances, with a risk of death.
- Its effects on the brain are unknown. Stimulation to the head or on each side of the skull must be avoided.
- Stimulation on the sides of the neck or on the carotid artery can cause serious adverse effects on your blood pressure or your heart rhythm.
- Stimulation on the front of the neck can cause severe muscle spasms that can block your airways and cause breathing difficulties.
- Do not position the electrodes on broken or injured skin, or which is dirty or unhealthy. Skin with irritation, sores or other lesions can lead to the injection of too much current on the area, which can cause burns.
- Do not place the electrodes near cancerous lesions because this may have a negative impact on these injuries.
- Do not place the electrodes on skin areas whose sensations are not normal. You may burn yourself due to a lack of feeling of the high intensity of the current.
- Do not put the electrodes on areas that are swollen, red, infected or inflamed or on skin rashes (eg. phlebitis, thrombophlebitis and varicose veins). Stimulation should not be performed on areas of thrombosis or thrombophlebitis because it can promote the circulation and lead to a greater risk of embolism.
- Do not put the electrodes on redness or open wounds. Open wounds may lead to applying too much current on the zone, causing burns. They can also favour the penetration of substances from the electrode into the skin.
- Do not place the electrodes on the inside of body cavities, such as in the mouth. Indeed, this device is only designed for external application.
- Do not make sudden movements during a session. This could cause a dysfunction of the device.

Do not use the device in the following conditions :

- Do not use the device if you are connected to high frequency surgical equipment. This could lead to burns on the skin under the electrodes and damage the device.
- Do not use the electro-stimulator if you are monitored by a doctor and you have not consulted him before using it.
- In the case of internal bleeding due to impacts or injury, do not use the device.
- To contract a muscle, do not use the electrical muscle stimulation in case of risk of muscle contraction that can disrupt the healing process. If the tendon or the muscle is torn, a muscle contraction can aggravate the wear, like a voluntary contraction. After recent surgery, after an acute trauma or fracture, this situation can also happen. In case of occurrence of a tendonitis, a muscle contraction can also aggravate the symptoms.
- Do not use the device while driving, operating the machines or any other activity during which the electrical stimulation may lead to a risk of injury.
- Do not use the appliance if you are subject to falling asleep during the session, as this may cause you to feel pain too late. If using at the time of goign to bed, set the timer so that the device does not turn itself off automatically.

- Never use **MyTens** in contact with water (in the bathroom, in the shower or in the pool, etc.) because this increases the risk of an electric shock and skin burns.

2. TRANSCUTANEOUS NEUROSTIMULATION INFORMATION

MyTens is a device designed for alleviating adults with acute post surgical pain or low to moderate chronic pain. It works on the principle of transcutaneous electrical neurostimulation (TENS), which enables you to relieve pain and soothe muscle tension. This device also helps to promote venous return and to strengthen the muscle mass. The device has been specially designed to be used at home. Transcutaneous electrical neurostimulation (TENS) is a non-analgesic drug therapy and without secondary effects, used by doctors and physical therapists for more than 30 years. The electrodes are placed on different areas of the body and transmit a low electrical current to the nerves. This stimulation brings the body to produce and disseminate endogenous analgesic-painkilling substances (endorphins, enkephaline) whose function is to anesthetize the pain. In case of pain or muscle tension, mini-electrical impulses in the nervous tissue may block the transmission of pain signals to the nervous system and trigger the release of endorphins. You can choose your programme among 11 preset TENS protocols on the **MyTens** application. For the strengthening of your muscle mass, you can select a programme among 8 protocols of electrical muscle stimulation (EMS). **MyTens** sends, via the electrodes placed on the area to stimulate, electrical pulses causing a muscle contraction which strengthens the muscle. It is controlled from the **MyTens** application. A diagram of the human body with treatment zones helps you to position the electrodes correctly (see paragraph 10). The **MyTens** application is available on the BewellConnect® platform of medical interpretation and custom follow-up downloadable free on the Store app or Google Play.

Counter indications

Do not use the device :

- In case of heart disease
- If you have epilepsy
- If you wear a pacemaker, a defibrillator or any other implanted or metal electronic device (for example a drug administration system). In these conditions, use could cause an electric shock, interference, burns or even death.
- If you are pregnant, during the first quarter. The effects of TENS on the development of the foetus are not yet known. During pregnancy, do not use the appliance on the uterus or abdomen to avoid triggering of contractions. Always consult your doctor or midwife if you are pregnant and intend to use the device.
- If you suffer from venous thrombosis or blood pressure, or thrombophlebitis. Using the electrical muscle stimulation (EMS) device could in this case cause a blood clot.
- In case of cognitive impairment
- If your pains have not been diagnosed, with the exception of the positive opinion of your doctor to use this device.

Side effects

- During the EMS session or shortly after, redness of the skin may occur at the level of the electrodes (around or below). They usually disappear within two hours of the session. In case of

persistence of the redness for more than 24h, please consult your doctor.

- An electrical muscle stimulation (EMS) session can cause muscle pain for some time.

Reactions

In the event of adverse effects in use, stop using this device and see your doctor. These side effects may occur in the form of:

- Nausea or fainting
- Burns at the location of the electrodes
- Painful sensations, headache
- Irritation of the skin in the place of the electrodes, even if the gel present on the electrodes does not present any known risk of allergic reactions to this day.

3. CHARACTERISTICS

- Pain relief and healing of muscle tension (TENS : 11 preset programmes)
- Strengthening of the muscular mass (EMS : 8 standard programmes)
- Stimulation of the venous return (TENS)
- Easy to use
- Data transfer onto a phone or a Bluetooth 4.0 tablet

4. USE

4.1 DESCRIPTION

See *schema A* page 2

- | | | |
|----------------------------------|--------------------|-----------------------|
| 1 On / off button | 2 Status indicator | 3 Removable connector |
| 4 Micro-USB socket of the device | 5 Electrodes | 6 Micro-USB cable |

4.2. PRECAUTIONS BEFORE USE


- Use the transcutaneous electrical neurostimulation programmes in bursts (TENS Burst) at rest. Indeed, they can cause a muscle contraction in the treatment area.
- In case of recent surgery, consult a doctor before using the device. The use of this device may affect the healing process.
- It is necessary to complete the treatment in the **MyTens** application or press the on/off button of the device before removing the electrodes. If you do not complete the treatment and remove the electrodes or the device, and if you touch the connectors, it is possible to feel pain in the fingers. This phenomenon, however, is not dangerous.
- Before any use of electrodes, check their expiration date on the packaging. If the date has expired, do not use them.
- Use new electrodes when the electrodes are damaged, dirty, less adhesive or if you begin to feel discomfort during stimulation (uncomfortable pinching, tingling).
- This product must only be used with the adapter, the cable and accessories recommended by the manufacturer.

- Before use, check that the device and the electrodes are not damaged. In case of any damage, do not use the device or one of the electrodes.
- Usable indoors and outdoors, the device is not designed to withstand all weather conditions.
- The device is not waterproof. Do not wet it and never use it in a damp environment.
- Do not use plaster or tape to attach the electrodes to the skin.
- For storing the electrodes, please follow the instructions in section 6.1 "Storage of electrodes".
- In case of abnormal operation of the device in accordance with this user manual, discontinue use and contact the after-sales department. Consult section 7 "Troubleshooting Guide" for information concerning potential malfunctions of the device.


4.3 PRESENTATION OF THE PROGRAMMES

4.3.1. FORMS OF WAVES


CONT : Form of continuous wave


HANS : Form of alternating wave

BURST : Wave form deviation

4.3.2. TENS PROGRAMMES

Traditional transcutaneous neurostimulation programme :

A common transcutaneous neurostimulation programme is performed with a high frequency stimulation. This stimulation is intended to block the transmission of the pain signal to the nervous system. Generally, people use traditional transcutaneous neurostimulation for a session of about 30 minutes. It is possible to use these programmes for a variable period and frequency according to what you want. It is normal to feel the effects of this programme more intensely during the stimulation and to experience its decay after the end of the session. To gain in efficiency, it is necessary to adjust the intensity of the programme, avoiding any painful or unpleasant sensation.

Burst transcutaneous neurostimulation programme :

A burst programme (Burst TENS) allows high frequency stimulation in the form of series of impulses. It enables the release of endorphins. It is possible to use the transcutaneous neurostimulation in bursts for 30 minutes. The frequency can be up to several times per day, 3 times for example. For a sensation of strong stimulation, it is advisable to adjust the intensity of the burst programme. An unpleasant feeling may then appear. You can also adjust the intensity of the programme at a level leading to perceptible muscle contractions in the stimulation area. These contractions are not dangerous. However, if you experience discomfort, it is advisable to decrease the intensity of the programme, change or move the electrodes. The burst transcutaneous electrical neurostimulation (Burst TENS) can trigger a muscle contraction in the area of stimulation. We therefore ask you to perform this programme at rest.

Summary Table of the TENS Programmes

| TENS | | | | Impulses | | | |
|----------------|--|----------------|------------------|------------------|-------------|-----------|---|
| Treatment site | Name | Treatment time | Treatment phases | Frequency (time) | Width in µS | Wave form | Description |
| Shoulder | Shoulder pain and wrinkled shoulder froissée | 30mn | Phase 1 | 2Hz (10s) | 250 | CONT | Continuous output |
| | | | Phase 2 | 4Hz (8s) | 250 | | |
| | | | Phase 3 | 6Hz (6s) | 250 | | |
| Lumbar areas | Pain in the lower back | 30mn | Phase 1 | 2Hz (10s) | 250 | CONT | Continuous output |
| | | | Phase 2 | 4Hz (8s) | 250 | | |
| | | | Phase 3 | 6Hz (6s) | 250 | | |
| | | | Phase 4 | 2Hz (10s) | 250 | | |
| | | | Phase 5 | 4Hz (8s) | 250 | | |
| | | | Phase 6 | 6Hz (6s) | 250 | | |
| Elbow | Pain in the elbow | 30mn | Phase 1 | 2Hz (3s) | 150 | HANS | Alternating output |
| | | | Phase 2 | 100Hz (3s) | 200 | | |
| Forearm | Pain in the forearm | 30mn | Phase 1 | 2Hz (3s) | 200 | HANS | Sortie alternative |
| | | | Phase 2 | 100Hz (3s) | 150 | HANS | |
| | | | Phase 3 | 3Hz | 250 | EMS | Up : 1,5s Down : 3s 10s cycle with 2/5 strength |
| Wrist | Pain in the wrist and carpal tunnel syndrome | 30mn | Phase 1 | 2Hz | 250 | CONT | Continuous output |
| Hip | Pain in the hip | 30mn | Phase 1 | 6Hz (30s) | 250 | CONT | Continuous output |
| | | | Phase 2 | 6Hz (30s) | 250 | | |
| | | | Phase 3 | 8Hz (20s) | 250 | | |
| | | | Phase 4 | 8Hz (20s) | 250 | | |
| | | | Phase 5 | 10Hz (20s) | 250 | | |
| | | | Phase 6 | 10Hz (20s) | 250 | | |
| Thigh | Sciatic pain | 30mn | Phase 1 | 100Hz | 200 | EMS | Up : 0,5s Con : 7s Down : 0,5s Rest : 7s |
| Knee | Knee pain | 30mn | Phase 1 | 100Hz | 150 | BURST | On : 0,25s Off : 0,25s |

| | | | | | | | |
|-----------|---|------|---------|-----------|-----|-------|---------------------------|
| Lower leg | Pain in the lower leg | 30mn | Phase 1 | 100Hz | 150 | BURST | On : 0,25s Off : 0,25s |
| Ankle | Pain in the ankle and Achilles tendon problem | 30mn | Phase 1 | 100Hz | 200 | CONT | Continuous output |
| Foot | Pain in the foot | 30mn | Phase 1 | 40Hz (5s) | 250 | CONT | Continuous output |
| | | | Phase 2 | 6Hz (10s) | 250 | | |
| | | | Phase 3 | 50Hz (5s) | 250 | | |

4.3.3. ELECTRICAL MUSCLE STIMULATION (EMS) PROGRAMMES

These programmes, commonly called “neuromuscular electrical stimulation” (NMES) or “electromyostimulation”, emit impulses to cause a muscle contraction. They are especially used for strengthening the muscle mass. Only perform the electrical muscle stimulation programme on intact and healthy muscles. Because they involve muscle contraction, they need to be performed at rest (lying or sitting).

When using for the first time, the electrical stimulation can lead to a atypical and strange sensation. It is advisable to begin the programme on low intensity and change to a higher one to become used to it. The stimulation must remain comfortable and muscle contractions must not be harsh or painful. If your programme is part of intensive physical exercise or after effort, you are advised to choose a low intensity to limit muscular fatigue.

Summary Table of the EMS Programmes

| EMS | | | Impulses | | | | |
|----------------|-----------------------------|----------------|------------------|------------------|-------------|-----------|--|
| Treatment site | Name | Treatment time | Treatment phases | Frequency (time) | Width in µS | Wave form | Description |
| Shoulder | Resistance of the trapezoid | 28mn | Phase 1 | 5Hz (5mn) | 280 | EMS | Preparation phase |
| | | | Phase 2 | 55Hz (10s) | 280 | | Intensity : 50% Up : 1,5s Down : 0,75s |
| | | | Phase 3 | 6Hz (8s) | 280 | | Intensity: 25% |
| | | | Phase 4 | 3Hz (10mn) | 280 | | Phase of relaxation Intensity : 40% |

| | | | | | | | |
|--------------|---------------------------------|------|---------|------------|-----|-----|---|
| Abdominal | Resistance of abdominal muscles | 22mn | Phase 1 | 30Hz (2mn) | 200 | EMS | Up : 10s Con : 5s Down : 10s Rest : 5s |
| | | | Phase 2 | 45Hz (9mn) | 200 | | Up : 5s Con : 5s Down : 5s Rest : 5s |
| | | | Phase 3 | 60Hz (9mn) | 200 | | Up : 6s Con : 8s Down : 6s Rest : 5s |
| | | | Phase 4 | 30Hz (2mn) | 200 | | Up : 10s Con : 5s Down : 10s Rest : 5s |
| Lumbar areas | Resistance of the lower back | 28mn | Phase 1 | 5Hz (5mn) | 300 | EMS | Preparation phase |
| | | | Phase 2 | 55Hz (10s) | 300 | | Intensity : 50% Up : 1,5s Down : 0,75s |
| | | | Phase 3 | 6Hz (8s) | 300 | | Intensity: 25% |
| | | | Phase 4 | 3Hz (10mn) | 300 | | Phase of relaxation Intensity : 40% |

| | | | | | | | |
|---------|--|------|------------|-------------|---------------------------------------|-----|---|
| Forearm | Front-arm : Extender of the wrist and fingers | 28mn | Phase 1 | 5Hz (5mn) | 300 | EMS | Preparation phase |
| | | | Phase 2 | 55Hz (10s) | 200 | | Intensity: 50% Up: 1,5s Down: 0,75s |
| | Phase 3 | | 6Hz (8s) | 200 | Intensity: 25% | | |
| | Phase 4 | | 3Hz (10mn) | 200 | Phase of relaxation Intensity: 40% | | |
| Hip | Resistance of the hip muscle | 32mn | Phase 1 | 5Hz (5mn) | 300 | EMS | Preparation phase |
| | | | Phase 2 | 75Hz (6,3s) | 300 | | Intensity: 50% Up: 1,5s Down: 0,75s |
| | | | Phase 3 | 4Hz (8s) | 300 | | Intensity: 25% |
| | | | Phase 4 | 3Hz (10mn) | 300 | | Phase of relaxation Intensity: 40% |
| Thigh | Resistance of the hamstring muscle | 28mn | Phase 1 | 5Hz (5mn) | 370 | EMS | Preparation phase |
| | | | Phase 2 | 55Hz (10s) | 370 | | Intensity: 50% Up: 1.5s Down: 0.75s |
| | | | Phase 3 | 6Hz (8s) | 370 | | Intensity: 25% |
| | | | Phase 4 | 3Hz (10mn) | 370 | | Phase of relaxation Intensity: 40% |

| | | | | | | | |
|-----------|--------------------------------|------|---------|------------|-----|-----|---|
| Lower leg | Resistance of the calf muscles | 28mn | Phase 1 | 5Hz (5mn) | 370 | EMS | Preparation phase |
| | | | Phase 2 | 55Hz (10s) | 370 | | Intensity: 50% Up: 1.5s Down: 0.75s |
| | | | Phase 3 | 6Hz (8s) | 370 | | Intensity: 25% |
| | | | Phase 4 | 3Hz (10mn) | 370 | | Phase of relaxation Intensity: 40% |
| Foot | Training the foot | 30mn | Phase 1 | 50Hz | 200 | EMS | Up: 2s Con: 8s Down: 2s Rest: 4s |

4.4. INDICATIONS BEFORE USE

Charging the device

Firstly charge the device before use. In normal times (ambient temperature), it takes 5 hours to charge the device. A full charge allows a battery life of about 8 hours continuously when functioning in normal conditions. Please charge the device fully before beginning a session. If the charge is not complete at the beginning of a programme, the battery may become depleted before the end of the session. You cannot use the device when it is on charge.

See Schema B page 2

1. Disconnect the cable used to connect the two units of the micro-USB of the unit with the battery light .
2. Insert the micro-USB cable into the micro-USB socket for the respective unit .
3. Insert the USB plug into a computer to charge the device.
4. When the battery light goes into continuous orange (which means that the battery is being charged) to a continuous green (the battery is full), remove the micro-USB plug from the micro-USB socket of the unit with the battery indicator (see diagram 4). Use the USB plug only for the use described above (to charge and connect to the cable connecting the two units).
5. Reconnect the cable connecting the two units : Insert the connector of this cable into the micro-USB socket for the unit with the battery light .

4.5. DOWNLOADING THE APPLICATION

Download the BewellConnect® application to your phone or tablet:

- by the App store or Google Play

- or by flashing the QR Code on the side of the box.

Then click the **MyTens** icon.

4.6. STARTING AND USING MYTENS

Preparation of the treatment areas and installation of myTens

1. Before placing the electrodes, check that the skin does not show injury or redness and is healthy.
- Transcutaneous electrical neurostimulation (TENS) programmes : Do not position the electrodes on irritated skin or skin with open wounds, rashes or cancerous lesions
- Electrical muscle stimulation (EMS) programmes : Position the electrodes only on intact and healthy muscles.
2. Using a damp cloth, clean the skin and then dry it. To receive the electrodes, it must be clean and dry, without cream or lotion. Any oil, dust or other element that may be on the skin could prevent the electrodes from adhering.
3. Remove the electrodes from the bag, ensure they are not damaged or worn. If the expiry limit is exceeded, do not use these electrodes.
4. Attach the electrodes to the device using the pressure buttons provided for this purpose.
5. Put the first electrode on the treatment area chosen. See the section "Arranging the electrodes" and the electrode positioning diagrams at the end of this manual (see paragraph 10).
6. Put the second electrode on the area of the body to treat (Voir paragraphe 10). Avoid any contact between the electrodes and not place one on the other.
7. To apply the second pair of electrodes, proceed in the same manner (steps 1 to 6). The second pair must be placed symmetrically (e.g. for treatment of the left thigh with one pair of electrodes, the second pair must be applied symmetrically on the right thigh). The intensity of the electrical impulses will be identical in both pairs, unless otherwise selected.

Remarks

- It is possible to perform a session dressed. Simply place the electrodes and the devices fixed under clothing.
- You can keep a device on the body, even between two sessions.
- It is possible that you may need external assistance for placing the electrodes in certain areas.

Locating the electrodes

Warning


As indicated in "Do not treat the the following parts or areas" :

Do not place the electrodes on the head, carotid, on any side of the neck, at the front of the neck or on the torso. Do not position the electrodes on swollen skin or skin with open wounds, rashes or cancerous lesions or on irritated skin.

The effectiveness of your sessions will depend in part on how you place your electrodes on the skin. You are advised to place the electrodes at the level of the area to be treated or in the vicinity of this area. Before installing the electrodes in a place, ensure you have good feeling in the area. If in doubt about the ideal place, do not hesitate to move the electrodes to find the best feeling. In case of proven effectiveness on a specific area, do not hesitate to note the place for using there again. If you suffer discomfort (unpleasant contractions during burst electrical neurostimulation programmes), it is advisable to move the electrodes. If during electrical muscle stimulation programmes (intended especially to create a contraction of the muscle) no contraction occurs, this may mean an intensity too low or bad installation of the electrodes.

Starting a session through the use of the MyTens application


1. Ensure that the Bluetooth (4.0) connection is activated on your smartphone or tablet. Open the **MyTens** application .
2. Press the on/off button so that the appliance is put on standby mode (Continuous green light : Standby mode, ready for the session). The App detects the **MyTens** which is switched on. In case of inactivity for more than 5 min, the device switches off. Press the on/off button again to begin a session.
3. Choose between a TENS session (transcutaneous electrical neurostimulation) or EMS session (neuromuscular electrical stimulation).

4. In the case of past use of **MyTens**, you can start again the last session or choose another one in the history.


History Button


Select the part to be treated

In the case of past use of MyTens or a new session, it is essential to select a treatment area and a program.


Estimate the intensity of pain (in TENS mode)

On a scale of 1 to 10, evaluate the intensity of the pain. For no pain press 0, for extremely intense pain press 10.

Select and start a session.


1. Starting a session : Click on Play Button.

2. Stopping a session : Click on the Stop Button.

3. Estimate your satisfaction level (the satisfaction level is available in TENS mode only). The session is then recorded in the history.

4. To turn off the device, press the “On/Off” button.

5. Remove the electrodes from the skin. Detach from the edges.

6. Remove the device from the electrodes. If you forget the device for more than 5 minutes, it turns off.

If the device is frequently used during the day, it is advisable to charge it before starting a new session.

5. STORAGE AND MAINTENANCE

5.1. STORAGE OF THE ELECTRODES

Put the electrodes back into their protective bag. To prevent buildup of dust on the electrodes, close the bag carefully. Prefer to store electrodes in a place where the temperatures is 5 to 27°C. Avoid extreme heat and exposure to direct sunlight. Do not keep the electrodes in the freezer or the refrigerator.

5.2. CLEANING

Warning : Never immerse the product in water and do not rinse it. Never immerse the electrodes in water. Do not use cleaning agents other than those mentioned below. They could cause serious damage to the equipment.

1. Using a damp cloth and mild detergent (eg: dishwashing liquid), clean the device. It is also possible to use isopropyl alcohol at 70° (IPA).

2. If dirt gets on the electrodes, put a drop of water on your finger and gently remove the dust from the surface. Do not use soap or alcohol to clean the electrodes.

6. TROUBLESHOOTING GUIDE

| Problem | Possible cause | Solution |
|---|---|--|
| The battery indicator on the device is displayed in orange continuously | Discharged battery | Recharge the device |
| Orange battery light and then the device switches off. | Discharged battery | Recharge the device |
| Treatment which seems to be different from or less pleasant than the previous one | <ol style="list-style-type: none"> 1. Electrodes not placed in the correct location 2. Intensity too low or too high 3. Expiry date on the electrodes exceeded 4. Worn electrodes 5. Electrodes too dirty 6. You did not correctly apply the electrodes to the skin | <ol style="list-style-type: none"> 1. Switch off the device and remove the electrodes. Remove the electrodes and move them slightly. 2. In the App, modify the intensity level. 3. Remove the electrodes and replace them with new ones. 4. Replace the electrodes. 5. Clean the electrodes (paragraph 6.2 “Cleaning”). If the problem occurs again, replace the electrodes. 6. Switch off the device and remove the electrodes. Make sure the electrodes are properly placed on the skin. |
| Unpleasant feeling when touching the connectors of the device or the electrodes | The session has not been stopped in the application or the on/off button has not been pressed before removing the device | Always end the session in the application or remember to press the on/off button before removing the device. |

7. AVOIDING MALFUNCTIONS

- Do not perform stimulations in the vicinity of electronic surveillance devices (such as cardiac monitors or electrocardiogram alarms). They might not work properly when using the electrical stimulation device at the same time.

- Do not operate the device at less than 1 metre from shortwave or microwave medical equipment. Close proximity with this device may cause a degradation of performance of the device.

- Portable RF communications equipment (including devices like antenna cables and external antennas) should not be used at less than 30 cm from the least part of the device, including cables mentioned by the manufacturer. This could lead to deterioration of performance of the device.






-The use of accessories, transducers and cables other than those mentioned or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and increase the risk of malfunction






8. TECHNICAL CHARACTERISTICS

| | |
|---------------------------------|---|
| Units | 4 devices (49,9mmx55,8mm) |
| Programmes | 19 programmes (11 TENS + 8 EMS programmes) |
| Battery | 3.7V/500mAh Li-Ion |
| Output voltage | max 60mA ±10% (1000Ω load) |
| Conditions of use | 5°C to 40°C with a relative humidity of 15% to 93%, atmospheric pressure of 700 to 1060 hPa |
| Storage conditions | -10°C to 50°C with a relative humidity of 10% to 90%, atmospheric pressure of 700 to 1060 hPa |
| Dimensions | 440x54x13mm cable included |
| Weight | 55g (±20%) by unit |
| IP classification system | IP22 |
| Wave form | Biphasic square wave |
| Impulse width | 150us~370us (±10%) |
| Heart rhythm | 2Hz~100Hz (±10%) |
| Treatment time | 1-90 minutes (±10%) given by the application |
| Adapter | 1pc (Input : 100-240 V AC 50/60Hz Output : 5 V DC 300mA) |







9. PLACING THE ELECTRODES

9.1. POSITION OF THE ELECTRODES-TENS PROGRAMMES

| | | |
|--------------|--|--|
| Shoulder | Shoulder pain and wrinkled shoulder |  |
| Lumbar areas | Pain in the lower back |  |
| Elbow | Pain in the elbow |  |
| Fore-arm | Pain in the forearm |  |
| Wrist | Pain in the wrist and carpal tunnel syndrome |  |
| Hip | Pain in the hip |  |

| | | |
|-----------|---|---|
| Thigh | Sciatic pain |  |
| Knee | Knee pain |  |
| Lower leg | Pain in the lower leg |  |
| Ankle | Pain in the ankle and Achilles tendon problem |  |
| Foot | Pain in the foot |  |

9.2. POSITION OF THE ELECTRODES-EMS PROGRAMMES

| | | |
|--------------|--|--|
| Shoulder | Resistance of the trapezoid |  |
| Abdominal | Resistance of abdominal muscles |  |
| Lumbar areas | Resistance of the lower back |  |
| Fore-arm | Extender of the wrist and fingers Bender of the wrist and fingers |  |
| Hip | Resistance of the hip muscle |  |
| Thigh | Resistance of the hamstring muscle |  |

| | | |
|-----------|--------------------------------|--|
| Lower leg | Resistance of the calf muscles | |
| Foot | Training the foot | |

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EN: Visiomed* will repair or replace this product free of charge in the case of defective parts or manufacturing defects, in accordance with the conditions mentioned below as follows:

DURATION: 24 MONTHS RETURN TO WORKSHOP

LIMITS AND EXCLUSIONS: This guarantee concerns only the original final purchaser. A purchase invoice, or another proof of purchase, with this guarantee card will be required to obtain an after-sales service, in accordance with this guarantee. This guarantee card will not be extended to another person only the original final purchaser. This guarantee becomes void if the serial numbers on the product are modified, replaced, illegible, absent, or if repair has been carried out by a service not approved, including the user.

This guarantee covers only the defects of the material or parts, occurring during normal use of the product. It does not cover the damage caused during the transport of the apparatus, caused due to repairs being carried out by the distributor, by any modifications undertaken, any connection of equipment not approved by Visiomed*, or causes contrary to those written in the user manual or notice. Moreover, the present guarantee does not cover damage due to falls, bad handling, bad installations, damage by fire, floods, lightning, or any other natural disaster. This guarantee does not cover the packing of the material, the accessories, the defects caused by commercial exposure of the product, show room, sale space, demonstration etc... Normal maintenance, cleaning and the replacement of parts where wear is normal, are not covered by the terms of this guarantee. Visiomed* and its representatives and agents will not in any case be held responsible for any damage and consecutive damages due to the mishandling of this product. This guarantee is the only valid one at Visiomed*, any other guarantee (commercial guarantee) except this one will not be taken into account.

IMPORTANT: During the guarantee period if you are dissatisfied with the repairs of this product, please contact the Visiomed* customer service.

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EN: *Disposal of Waste Equipment by Users in Private Households in the European Union.*

The symbol on the product or on its packaging indicates that this product must not be disposed of with your other household waste. Instead, it is your responsibility to dispose of your waste equipment by handing it over to a designated collection point for the recycling of waste electrical and electronic equipment. The separate collection and recycling of your waste equipment at the time of disposal will help to conserve natural resources and ensure that it is recycled in a manner that protects human health and the environment. For more information about where you can drop off your waste equipment for recycling, please contact your local city office or your household waste disposal service.

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EN: *BewellConnect® is a brand registered by Visiomed Group SA. The constructor keeps the right to modify the photographs and characteristics without previous notice. Registered design. The logos and trademarks are the property of their respective owners.*

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© 07/2015 BW-TSX