ThermoFlash®

EC REP
VISIOMED GROUP SA
112, Avenue Kléber
75116 PARIS - France

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
ThermoFlash®

USER MANUAL

MicroSecond Flash™
TECHNOLOGY
NEW GENERATION

THERMOFLASH® LX-361T is clinically tested and approved ASTM 1965-1998 (2009).

CONTENTS

I. INTRODUCTION.....	19
II. PRECAUTIONS TO TAKE BEFORE USE.....	19
III. USE.....	19
DIFFERENT MEANS OF READING BODY TEMPERATURE.....	20
ADVANTAGES OF TAKING THE TEMPORAL ARTERY TEMPERATURE (TA).....	20
GUIDELINES TO FOLLOW BEFORE TAKING A TEMPERATURE.....	20
HOW TO TAKE A TEMPERATURE.....	21
CONSTRAINTS.....	21
IV. FUNCTIONS.....	21
V. DESCRIPTION OF THERMOFLASH® LX-361T.....	21
VI. STARTING UP YOUR THERMOFLASH® LX-361T.....	22
VII. CONFIGURATION & MENU FUNCTIONS.....	22
1. CHOOSING HOW TO TAKE THE TEMPERATURE.....	22
2. CHOOSING THE TEMPERATURE UNIT.....	22
3. THE MEM BUTTON.....	22
F1 – HIGH TEMPERATURE ALERT SETTING.....	22
F2 – RECALIBRATION.....	22
4. EXIT CONFIGURATION MODE.....	23
5. THE BACKLIGHT BUTTON.....	23
6. THE LANGUAGE BUTTON.....	23
7. DATA SAVING.....	23
8. CHANGING THE BATTERIES.....	23
VIII. PRECISION & TECHNICAL CHARACTERISTICS.....	23
LONGEVITY USE.....	24
IX. NOTES.....	25
X. ACCESSORIES.....	25
XI. TROUBLESHOOTING.....	25
GUARANTEE AT THE END OF THE MANUAL.....	

- Follow the maintenance instructions listed in this document.
- This thermometer may be used for personal use at home.
- Only use this thermometer for the prescribed purposes listed in this document.
- Use this thermometer in an ambient temperature between 10°C and 40°C.
- This thermometer must always be stored in a clean, dry place.
- Do not expose this thermometer to electric shocks.
- Do not expose this thermometer to extreme temperature conditions: T° > 50°C, T° < -20°C.
- Do not use this thermometer in an environment with higher than 85% relative humidity.
- The glass protecting the lens is the most fragile part of the thermometer.
- Do not touch the infra-red lens with your fingers.
- Clean the glass with a cotton bud dipped in 70% alcohol.

- Do not expose the thermometer to direct sunlight or to water.
- Do not use this thermometer outdoors.
- Never drop the thermometer.
- Do not attempt to repair this thermometer. In the event of a problem, contact your retailer.
- Please keep the thermometer off the children's touch. Prevent the product from children playing, swallowing and eating.
- The MEDICAL ELECTRICAL EQUIPMENT needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the ACCOMPANYING DOCUMENTS.
- Portable and mobile RF communications equipment can affect MEDICAL ELECTRICAL EQUIPMENT.
- The measuring result of this product is only for your reference. If you have any doubt, please measure the temperature in other methods.

IMPORTANT

Check the batteries are inserted the right way. Reversing the polarity may damage the thermometer and affect its warranty.

The manufacturer reserves the right to alter the specifications of the product without prior notification.

I. INTRODUCTION

The **THERMOFLASH® LX-361T** is a revolutionary infra-red medical thermometer, fitted with cutting edge MicroSecondFlash™ automatic calibration technology, developed by **Visiomed®**.

Easy to use, it reads your body temperature from a distance of 2 to 5 cm in less than a second. Fast and accurate, all you have to do is position the **ThermoFlash®** toward the forehead to instantly see the temperature. **ThermoFlash®** is especially effective for taking the temperature of babies and small children. It has been demonstrated that for new-born children, taking the temperature via the temporal artery is more accurate than taking it aurally, and less distressing than via the rectum*.

ThermoFlash® can be used at any age.

ThermoFlash® can also be used to check the temperature of the bath, food, or ambient temperature.

* *Greenes D, Fleisher G. Accuracy of a Non-invasive Temporal Artery Thermometer for Use in Infants. Arch Pediatr Adolesc Med 2001; 155:376.*

II. PRECAUTIONS TO TAKE BEFORE USE

The parameters for the **THERMOFLASH® LX-361T** are set at the factory. It is not necessary to calibrate the thermometer during start-up.

In order to guarantee a reliable and stable result, we advise when moving to an environment of a very different temperature, to leave your **THERMOFLASH® LX-361T** at this ambient temperature for 15-20 minutes before use. Before taking the next temperature, wait for the screen to switch off. You cannot take a temperature while the screen is lit up.

III. USE

Every solid, liquid or gaseous body emits energy by radiation. The intensity of this energy depends on the temperature of the object. The infra-red thermometer **THERMOFLASH® LX-361T** can in this way read a person's temperature by the quantity of energy he or she emits. This is possible thanks to the external temperature probe which constantly analyzes and records the ambient temperature. This means that as soon as the operator sets the thermometer close to the body and activates the heat sensor, the reading takes place instantly, detecting the infra-red heat emitted by the arterial blood flow. The reading

is thus obtained with no interference whatsoever from ambient heat. **ThermoFlash®**'s high quality is connected to this unique technology developed by **Visioned®**.

DIFFERENT MEANS OF READING BODY TEMPERATURE

- **Core temperature:** Core temperature gives the most accurate reading. It relies on measuring the temperature of the pulmonary artery via a catheter fitted with a thermal probe, giving the temperature in situ. The same principle is used for probes measuring the temperature in the esophagus. However, such invasive temperature measurement method requires specific equipment and expertise.

- **Rectal temperature:** Rectal temperature adjusts slowly in comparison to the evolution of body temperature: it has been demonstrated that rectal temperature remains high even after the patient's internal temperature has started to drop and vice-versa. Also, using this method, rectal perforations have been reported, and, without the appropriate sterilization techniques, rectal temperature readings may propagate the spread of contaminants often found in faeces.

- **Oral temperature:** Oral temperature readings are easily influenced by recent consumption of food or drink and by breathing through the mouth. To measure the oral temperature, the mouth must remain closed and the tongue lowered for three or four minutes. This is not easy to achieve with young children.

- **Axillary temperature:** It may be easy to take a reading of the axillary temperature, but it has been demonstrated that it gives an inexact representation of the child's internal temperature. To take this kind of reading the thermometer must be wedged tightly over the axillary artery. Despite its relatively low sensitivity in detecting fevers, the axillary temperature reading is recommended by the *American Academy of Pediatrics* as a screening test for a fever in a new-born child.

- **Aural temperature:** It is essential to take the aural temperature in the correct way to ensure an accurate reading. The thermometer probe must be against the warmest part of the inner ear. Failure to do so may result in an inaccurate reading.

ADVANTAGES OF TAKING THE TEMPORAL ARTERY TEMPERATURE (TA)

Infrared arterial temperature can be measured by using a device positioned on the forehead, near the temporal artery. It has been demonstrated that this relatively new method is more precise than taking the temperature aurally, and is less distressing than via the rectum, especially for new-born children.

The **THERMOFLASH® LX-361T** thermometer has been designed to give an immediate and no contact reading from the temporal artery temperature. This artery is close enough to the surface of the skin to be easily accessible, and gives an accurate reading since the blood flow is permanent and regular. This is connected to the heart by the carotid artery, which is directly connected to the aorta. It is part of the core arterial network. The efficiency, speed and comfort of taking a temperature reading from this area makes it ideal compared to other methods.

GUIDELINES TO FOLLOW BEFORE TAKING A TEMPERATURE

- In order to ensure an accurate and reliable reading, it is essential that each user is suitably informed and knows how to take a temperature with this thermometer.

- It is essential to remember that although procedures such as taking temperature may be simple, they should not be trivialized.

- Temperature should be taken in a neutral context. The patient must not be recovering from an intense physical activity before the measurement, and the ambient temperature must be moderate.

- Be aware that physiological variations should be taken into account when evaluating the results: the temperature increases by 0.5°C between 6am and 3pm. Women in general are 0.2°C warmer than men,

and their temperature varies according to the menstrual cycle. It increases by 0.5°C in the second phase of the cycle and at the beginning of pregnancy.

- While sitting, the temperature is 0.3 to 0.4°C lower compared to while standing.

HOW TO TAKE A TEMPERATURE

SEE SCHEMA **A** PAGE 2.

Aim at the FOREHEAD, over the right temporal region, from a distance of about 2-5cm (0.8-2 in), press the thermometer's measurement button and the temperature is instantly displayed.

The reliability of the measurement cannot be guaranteed if the temperature is measured over another part of the body (e.g. arm, torso...).

CONSTRAINTS

In order to ensure a stable and reliable reading, before each measurement:

- Push back hair from the forehead.

- Wipe away any perspiration from the forehead.

- Avoid any drafts (e.g. from nasal specs, air conditioning...).

- Allow a **1 minute interval** between two measurements.

- Any time you move the **ThermoFlash®** to a place with a significantly different ambient temperature, **leave it to acclimatize at this temperature for at least 15 minutes before use.**

IV. FUNCTIONS

1. Specially designed for a body temperature reading without interference from ambient heat.

2. Reliable and stable results thanks to the HEIMANN infra-red detection system.

3. Vocal indication of measurements and color indicators to quickly display and interpret the diagnosis of the temperature obtained.

4. An alert if the temperature is too high.

5. A memory unit to save the last 32 readings.

6. Backlit digital LCD screen.

7. Choice of reading in Celsius or Fahrenheit.

8. Automatic shut-down (to save energy).

9. Small, practical and easy to use.

Secondary functions: The **THERMOFLASH® LX-361T** also lets you take the temperature of a bath, a bottle or the ambient temperature of a room (in **SURFACE TEMP** or **ROOM** mode).


V. DESCRIPTION OF THERMOFLASH® LX-361T

SEE SCHEMA **B** PAGE 2.

- | | | |
|-------------------------|--------------------------|----------------------------|
| 1 Infra-red lens | 4 Battery case | 6 LED for night use |
| 2 LCD screen | 5 Setting buttons | 7 Infra-red sensor |
| 3 Measure button | | |

VI. STARTING UP YOUR THERMOFLASH® LX-361T

- Insert the batteries.
- At the first use or after changing the batteries, wait around 10 minutes for the thermometer to warm up.
- Aim at the forehead (SEE SCHEMA **A** PAGE 2) at 2-5 cm (0.8-2 in), press the measure button and the temperature will be displayed immediately. The temperature can also be taken behind the ear on the jugular.
- Before taking a reading, make sure there is no hair or sweat on the forehead.
- Change the batteries when shows on the screen.

VII. CONFIGURATION & MENU FUNCTIONS

1. CHOOSING HOW TO TAKE THE TEMPERATURE

The THERMOFLASH® LX-361T is specially designed to take a human body temperature with the **BODY** function.

You can also use the **ThermoFlash®** to measure the temperature of the surface of an object, a food or a liquid by selecting the **SURFACE TEMP** mode. By selecting the **ROOM** mode, it will display the ambient temperature of the room you are in.

Temperature range for **BODY** mode: 32°C - 42.9°C (89.6°F - 109.2°F).

Temperature range for **SURFACE TEMP** mode: 0°C - 60°C (32°F - 140°F).

Temperature range for **ROOM** mode: 0°C - 50°C (32°F - 122°F).

Press the *Measure* button then the *MODE* button for one second, the screen will display **BODY**.

Press the *MODE* button again to switch to **SURFACE TEMP** or **ROOM** mode.

Note: The THERMOFLASH® LX-361T is pre-calibrated to BODY.

Important: Surface temperature is not the same as internal body temperature. Internal temperature must always be read with the **BODY** function.

Always select the **BODY** function for an internal temperature reading and the **SURFACE TEMP** or **ROOM** function for an external reading (e.g. bottle, bath, bedroom...).

2. CHOOSING THE TEMPERATURE UNIT

Press the *°C/°F* button for one second: the THERMOFLASH® LX-361T will display the temperatures in degrees Celsius (the screen will display °C).

Press the *°C/°F* again, and the THERMOFLASH® LX-361T will display the temperatures in degrees Fahrenheit (the screen will display °F).

3. THE MEM BUTTON

F1. High temperature alert setting

Press and hold the *MEM* button for 3 seconds, and the thermometer will beep and the screen will show: **F1**. The alert setting will flash on and off.

Press the *°C/°F* button to lower the setting by 0.1°C, or the button *MODE* button to increase it by 0.1°C.

To save the setting, press the *MEM* button.

Note: the high temperature alert is by default set to 38°C.

F2. Recalibration

To recalibrate the temperature readings of your THERMOFLASH® LX-361T, press and hold the *MEM* button

for 3 seconds. The thermometer will beep and display: **F1**. Press the *MEM* button again to display the **F2** menu. Press the *°C/°F* button to lower the setting by 0.1°C, or the *MODE* button to increase it by 0.1°C.

To save the setting, press the *MEM* button.

If you are unsure, we recommend that you leave the setting at +0.0°C (default setting).

Any time you move the **ThermoFlash®** to a place with a significantly different ambient temperature, leave it at this temperature for at least 15 minutes before use.

4. EXIT CONFIGURATION MODE

Press and hold the *MEM* button until the screen switches itself off.

5. THE BACKLIGHT BUTTON

Press button 4 to activate the LED and to help you while taking a temperature at night.

6. THE LANGUAGE BUTTON

Use the language button to select the language for the voice function (French/English/Spanish/OFF).

The name of the language selected will be spoken. By default, the THERMOFLASH® LX-361T is set to French.

7. DATA SAVING

To display the last temperature readings saved, press the *MEM* button for 5 to 7 seconds.

You will then see the last temperature reading. Scroll through the readings with the *MODE* button.

The number displayed in the interval between two readings is the reading number.

Press the *°C/°F* button to scroll back.

To exit the memory menu, press the measure button. The **ThermoFlash®** will switch off automatically after five seconds.

8. CHANGING THE BATTERIES

Display: When the LCD screen displays **BATTERY** , the battery has run out.

Operation: Open the cover and replace the old batteries, making sure to insert them the correct way round. Inserting them incorrectly may damage the thermometer and void the thermometer's warranty. Never use rechargeable batteries. Use single use batteries.

VIII. PRECISION & TECHNICAL CHARACTERISTICS

1. Normal conditions of use operating temperature: 10°C- 40°C (50° - 104°F). Normal conditions of storage and transportation temperature: -20°C - 50°C (-68° - 122°F). Humidity rate: 10 - 85%

2. Power: DC3V (2 batteries LR03).

3. Size with base: 172,5 x 60,35 x 60,35 mm (L x W x H).

4. Weight: 120g.

5. Display Resolution: 0.1°C (0.1°F).

6. Measuring range:

In **BODY** mode: 32°C - 42.9°C (89.6°F - 109.2°F).

In **SURFACE TEMP** mode: 0°C - 60°C (32°F - 140°F) and **ROOM** mode: 0°C - 50°C (32°F - 122°F).

7. Energy consumption: ≤ 50mW.

8. Accuracy: 36°C - 39°C : ± 0.2°C.

9. Measurement distance: 2 cm - 5 cm (0.8 in - 2 in).

10. Automatic standby: 5 sec.

This product complies with the electromagnetic compatibility requirements EN60601-1-2 and the safety standards 60601-1, as specified in the EEC directive 93/42/EEC.

LONGEVITY USE

The THERMOFLASH® LX-361T is guaranteed for 40.000 readings.

EMC DATA

Guidance and manufacturer's declaration - electromagnetic immunity - for equipment and system that are not life-supporting. The model LX-361T is intended for use in the electromagnetic environment specified below. The customer or the user of the model LX-361T should assure that it is used in such an environment.

Immunity test	EN 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF	3 Vrms	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the model LX-361T, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [\frac{3.5}{E1}] \sqrt{P}$ 80 MHz to 800 MHz $d = [\frac{7}{E1}] \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). ^a Fields strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
IEC 61000-4-6	150 kHz to 80 MHz		
Radiated RF	3 V/m	3 V/m	
IEC 61000-4-3	80 MHz to 2.5 GHz		

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.

a. Fields strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model LX-361T is used exceeds the applicable RF compliance level above. The model LX-361T should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary such as reorienting or relocating the model LX-361T.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Recommended separation distances between portable and mobile RF communications equipment and the equipment or system - for equipment and systems that are not life-supporting.

The model LX-361T is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the model LX-361T can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the model LX-361T as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output of transmitter	Separation distance according to frequency of transmitter m	
	80 MHz to 800 MHz $d = [\frac{3.5}{E1}] \sqrt{P}$	800 MHz to 2.5 GHz $d = [\frac{7}{E1}] \sqrt{P}$
0.01	0.12	0.23
0.1	0.38	0.73
1	1.2	2.3
10	3.8	7.3
100	12	23

For transmitters rated at a maximum output power not listed above the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagations is affected by absorption and reflection from structures, objects and people.

IX. NOTES

- The glass protecting the lens is the most important and most fragile part of the thermometer, so be careful with it.
- Clean the glass with a cotton cloth, dampened with water or 70% alcohol.
- Only use the recommended batteries, do not recharge non-rechargeable batteries, do not burn them.
- Remove the batteries if the thermometer is not to be used for a long period.
- Do not expose the thermometer to direct sunlight or to water.
- A shock may damage the thermometer.

X. ACCESSORIES

User manual, Warranty card, Base, Batteries (LR03).

XI. TROUBLESHOOTING

If you have one of the following problems while using your THERMOFLASH® LX-361T, please refer to this breakdown service guide to help resolve the problem. If the problem persists, please contact our customer service at +33 892 350 334.

- **The screen is showing temperatures higher than 95°C:** the temperature is in Fahrenheit. Change the measurement unit to Celsius by pressing the Program setting button to get F1 function key.
- **The screen is showing temperatures lower than 32°C (89.6°F):** to take the internal temperature of your body, choose the BODY function. If you are in SURFACE TEMP or ROOM mode, the 32°C temperature displayed is showing the external temperature that your body releases.
- **The screen is showing the message HI:** the temperature taken is higher than the range of the mode selected, either higher than 42.9°C (109.2°F) in BODY mode or 60°C (140°F) in SURFACE TEMP mode.
- **The screen is showing the message LO:** the temperature taken is lower than the range of the mode selected, either lower than 32°C (89.6°F) in BODY mode or 0°C (32°F) in SURFACE TEMP mode.

This message appears in different circumstances - please find below the most common:

Causes of a Lo Message	Solutions
Reading hampered by hair, sweat...	Make sure there is no obstruction.
Reading hampered by an air flow...	Make sure you are not in a draught; this may interfere with the infrared system.
The reading was taken too close, the ThermoFlash® didn't have enough time to reset.	Respect a delay of 15 seconds between each reading - 1 minute is advised.
Measuring distance is too far.	Respect a measuring distance between 2 to 5 cm, away from the forehead.

ThermoFlash®


NL

HANDLEIDING

MicroSecond Flash™
TECHNOLOGY
NEW GENERATION

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EN: Visiomed* will repair or replace this product free of charge in the case of defective parts or manufacturing defects, in accordance with the conditions mentioned below as follows:

DURATION: 24 MONTHS RETURN TO WORKSHOP

LIMITS AND EXCLUSIONS: This guarantee concerns only the original final purchaser. A purchase invoice, or another proof of purchase, with this guarantee card will be required to obtain an after-sales service, in accordance with this guarantee. This guarantee card will not be extended to another person only the original final purchaser. This guarantee becomes void if the serial numbers on the product are modified, replaced, illegible, absent, or if repair has been carried out by a service not approved, including the user. This guarantee covers only the defects of the material or parts, occurring during normal use of the product. It does not cover the damage caused during the transport of the apparatus, causes due to repairs being carried out by the distributor, by any modifications undertaken, any connection of equipment not approved by **Visiomed***, or causes contrary to those written in the user manual or notice. Moreover, the present guarantee does not cover damage due to falls, bad handling, bad installations, damage by fire, floods, lightning, or any other natural disaster. This guarantee does not cover the packing of the material, the accessories, the defects caused by commercial exposure of the product, show room, sale space, demonstration etc... Normal maintenance, cleaning and the replacement of parts where wear is normal, are not covered by the terms of this guarantee. **Visiomed*** and its representatives and agents will not in any case be held responsible for any damage and consecutive damages due to the mishandling of this product. This guarantee is the only valid one at **Visiomed***, any other guarantee (commercial guarantee) except this one will not be taken into account.

IMPORTANT: During the guarantee period if you are dissatisfied with the repairs of this product, please contact the **Visiomed*** customer service.

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EN: Disposal of Waste Equipment by Users in Private Households in the European Union. The symbol on the product or on its packaging indicates that this product must not be disposed of with your other household waste. Instead, it is your responsibility to dispose of your waste equipment by handing it over to a designated collection point for the recycling of waste electrical and electronic equipment. The separate collection and recycling of your waste equipment at the time of disposal will help to conserve natural resources and ensure that it is recycled in a manner that protects human health and the environment. For more information about where you can drop off your waste equipment for recycling, please contact your local city office or your household waste disposal service.

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© 01/2017 LX-361T
EN: ThermoFlash® is a brand registered by Visiomed Group SA. The constructor keeps the right to modify the photographs and characteristics without previous notice. Registered design. The logos and trademarks are the property of their respective owners.